

Care Antibacterial Foam Soap Cartridge 1 Ltr



Art No	CHFCS00106
System of use	Hand Soap Foam
Packing type	Disposable cartridge
Content	1 Ltr
Country of origin	UAE

Antibacterial Foam Soap is a pearlised medicated cream hand soap foam. It is a viscous yellow liquid with a chloroxylenol – odor.

Why Care Antibacterial Foam Soap Cartridge ?

- › It removes all right to medium dirt and soilage from hands, yet is harmless to the skin.
- › It is mild, non allergic and safe for children.
- › Suitable for areas with high hygienic standards, such as medical facilities, public washrooms and food preparation areas.
- › Hygienic: the pump is disposed each time the cartridge is replaced.
- › It doses 15 ml of foam (0.5 ml of liquid) each time ensuring 2000 doses from this cartridge.
- › The disposable unit design is suitable with Multiflex dispenser which doesn't require service, cleaning or maintenance.
- › Good for normal traffic.

How to use Care Antibacterial Foam Soap Cartridge ?

Wet your hand with warm water, apply the foam, lather well. Rub hands together rapidly all over the hands and wrists, between the fingers, under the fingers nails and the back of the hands. Rinse well under running water and ensure all traces of foam are removed. Dry the hands with a towel or a warm air dryer.

Intercare Washroom & Hygiene guarantees the quality of the product . Its application recommendations are based on the test and practical experience .

we accept no responsibility for damage caused through improper use of the product .

Contact

Intercare LTD
Head Office || Sharjah-UAE || P.O.BOX: 5628
info@intercare.ae || sales@intercare.ae
Tel: +97165345500